

FLEX-RUN™ JUNIOR



The Flex-Run Junior is a custom designed solution for children. Its vertical compliance and efficient energy return are must-haves for high active children.

Designed to be the most flexible of the running feet, the Flex-Run Junior offers increased cushioning and lower-frequency dynamic makes this foot ideal for Transfemoral and Transtibial users.



USER INFORMATION

Amputation Level:	Transfemoral and Transtibial
Impact Level:	High to Extreme
Maximum Patient Weight:	84Kg (185lbs)



FOOT INFORMATION


Categories:	1-6
Sizes:	One Size
Weight of Foot:	336g (11.9oz)
Build Height Junior:	200mm (7 7/8")
Adapter Options:	Male or Female Bolt on Pyramid

CATEGORY SELECTION CHART

WEIGHT KG	15-24	25-34	35-45	46-56	57-67	68-84
WEIGHT LBS	33-53	54-75	76-99	100-123	124-147	148-185
High to Extreme Impact Level	1	2	3	4	5	6

FLEX-RUN™ JUNIOR

CLEARANCE




KITS

FLEX-RUN KIT


PART#
 FSX007

CATEGORY (1-6)




COMPONENTS

PART#	DESCRIPTION
FWX40001	Attachment Plate*



FOOT ADAPTER SETS

PICTURE	PART#	DESCRIPTION	WEIGHT	BOLT SIZE
1.	FWX00030	Male Pyramid Foot Adapter Set*	75g	M10 x 40mm
2.	CM150161	Female Pyramid Foot Adapter Set*	100g	M10 x 35mm



*Adult componentry is used for Flex-Run Junior

FLEX-RUN™ JUNIOR

KITS

REPLACEMENT BOLT

PART#	SIZE	CAT.
CM130221	M10	ALL
CM130224	M10	ALL